<i>Doc no</i>	<i>Doc title</i>	 LABRIDA
No. BC-016	Instructions For Use (IFU)	

LABRIDA

LABRIDA™ BioClean™

DA

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LABRIDA™ BioClean™

EN

Instructions for use of LABRIDA™ BioClean™

Please read this information carefully before using the product.

REF LBC2013.0001

BC-016 version 1

<p><i>Doc no</i></p> <p>No. BC-016</p>	<p><i>Doc title</i></p> <p>Instructions For Use (IFU)</p>	 <p>LABRIDA</p>
--	--	---

Product description

The Labrida BioClean is a dental instrument with a working end made of fast-degrading chitosan attached to a stainless steel stem covered with polypropylene. The LABRIDA BioClean reduces treatment time and allows for gentle and effective mechanical cleaning of dental implants. This product can be used in standard angled dental handpieces with speeds up to 2000 rpms. The best result is achieved using oscillating dental handpieces at 600-1000 cycles per minute.

The Labrida BioClean was developed as a disposable implant cleaning instrument. Due to the risk of cross-contamination and to achieve optimal cleaning results, the instrument must be replaced after cleaning each implant. The instrument cannot be sterilised for reuse.

Properties/functions

The Labrida BioClean allows effective cleaning of implants without damaging the implant surface. The Labrida BioClean chitosan instrument is gentle on tissue and does not leave residue that may induce a foreign body inflammatory reaction. The polypropylene covering of the stem protects against damage or discoloration of the prosthetic supra construction.

Indications

The Labrida BioClean was designed for the mechanical cleaning of exposed dental implant surfaces during regular maintenance treatment and to prevent inflammation (mucositis) around dental implants.

Instructions for use

Please comply with the basic guidelines for treatment in a sterile environment.

Open the part of the packaging covered by metal foil and carefully lift out the instrument grasping the metal holder of the instrument with a tweezers. Fasten the instrument to a dental handpiece and ascertain that the instrument is fastened properly to the handle. Handle the brush carefully prior to wetting and usage, and avoid touching the blue brush part with fingers or other instruments because this can damage the chitosan filaments. Clean the implant using the Labrida BioClean by carefully applying the instrument to the surface to be cleaned and moving the rotating instrument slowly across the surface of the implant. Make certain that the working head of the instrument is always touching the surface of the implant. Do not use force. We recommend irrigation with sterile saline. The recommended rotation speed for this instrument is 600-1500 rpms. Please use an oscillating dental handpiece for best results, if it is available. **This instrument is not designed for rotation speeds above 2000 rpms. If the rpms are too high, the instrument's working head will be destroyed, the stem may break and the instrument can detach from the metal holder.** Make certain that the brush does not wedge itself between the implant and surrounding tissue; this can injure the patient and damage the implant, or the stem can break or detach from the metal holder. Rinse the surface of the implant thoroughly with sterile saline after finishing the Labrida BioClean procedure.

Contraindications

The Labrida BioClean is not designed for use:

REF LBC2013.0001

BC-016 version 1

<i>Doc no</i>	<i>Doc title</i>	
No. BC-016	Instructions For Use (IFU)	

- on patients that have reported a reaction to or incompatibility with chitosan
- on clinically mobile implants
- on new implants that are still not integrated into the jawbone
- near large nerves and blood vessels if there is a risk of structural damage.

Side effects

The Labrida BioClean is made of chitosan (the working end), surgical steel and polypropylene (the stem). Chitosan seldom causes side effects, but the Labrida BioClean must not be used on patients who are allergic to or have experienced a negative reaction or incompatibility to chitosan, steel or polypropylene. Careless use of the instrument can injure the patient and damage the implant or equipment. Please read the instructions in the user's manual carefully and do not misuse the instrument.

Warnings and precautions

- The Labrida BioClean is a disposable sterile instrument.
- The instrument must be discarded after use due to the risk of cross-contamination and infection.
- The instrument must not be re-sterilised.
- The instrument must not be used if the package is not sealed properly or if it is damaged.
- The instrument must not be used if it is damaged.
- The instrument must not be used if the expiration date has passed.
- The Labrida BioClean must be stored at room temperature in a clean, dry and dark place.

The risk of complications can be minimised by correct patient selection, using good operating routines and through the correct use of the Labrida BioClean (including compliance with the instructions in the paragraphs on contraindications, warnings and precautions listed above).

Additional information

Please visit our web site (www.labrida.com) for additional and updated information about the BioClean instrument.

Validity

The validity of all previous versions of the user's manual and instructions expires upon the publication of any new version.

Ordering information

Labrida BioClean

BioClean is manufactured by LABRIDA AS

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